


**ORDER GUIDE**

1. Model should be packed properly including 
2. Before Shipping order must be informed & scheduled at **+91 8929023119**
2. We request to allow us 3-5 working days for quality work

Dr. Name ..... Account ID .....

Address .....  
City/State/ZIP .....

Patient ID/Name ..... ☐ Male ☐ Female Age .....

Deliver by 5 p.m. on ..... ☐ Call before starting case

Enclosed with case: ☐ Impressions ☐ Bite ☐ Models ☐ Articulator ☐ Shade Tab ☐ D-Wax ☐ Pre-Op Models ☐ Photos

**CHOOSE PRODUCT**


- ☐ Dmax strong  
☐ Dmax esthetic

**PONTIC DESIGN**



**\*\*Standard unless specified otherwise**

**INCISAL SHAPE INSTRUCTIONS**



- ☐ Rounded ☐ Squared ☐ Pointed

**OCCLUSAL STAINING**

- ☐ Light\*\* ☐ Med ☐ Dark ☐ None

**PRESENT TOOTH OR STUMP SHADE**



**Indicate Shade Here**


**FINAL CERAMIC SHADE**



**Indicate Shade Here**

**Online – **

**NOTE:** Please send a study model on all work involving anterior teeth.



## TERMS AND CONDITION

All the terms & Conditions valid only in INDIA





**TERMS:** Cost of collection to be paid by customer in case any. All orders are payable 50% before start of work, and balance to be paid against delivery or within 30 days of invoice. If Due is not cleared within time frame additional 2% collection charges will be applicable on time of collection.

**NO-FAULT REMAKE POLICY:** [www.vaardental.com](http://www.vaardental.com) and its associate franchise is pleased to offer all remakes or adjustments at no additional charges within the warranty period and accompanied by the return of original appliance

**LIMITED WARRANTY/LIMITATION OF LIABILITY.** For warranty terms and conditions and limitation of liability, visit [vaardental.com/policies-and-warranties/](http://vaardental.com/policies-and-warranties/).

## PREPARATION GUIDELINES



### **vaardental** Strong

- A. 1.0 mm ideal reduction (0.5 mm minimum)
- B. Chamfer or shoulder margins preferred. Feather-edge OK
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

### **vaardental** esthetic

- A. 1.25 mm ideal reduction (0.7 mm minimum)
- B. Chamfer or modified shoulder margins preferred
- C. Axial walls must be convergent (avoid undercuts)
- D. Preparation should be cut in three planes
- E. To achieve optimal impression quality, gingival retraction is necessary for preparations with subgingival or equigingival margins

## SIGNATURE ANTERIOR SERIES

Choose the desired esthetic outcome for your patient

### SQUARE-TAPERED



Exemplify boldness and youthfulness.

### SOFT-SQUARED



Show energetic professionalism.

### OVOID



Convey charm and softness.

### TRAPEZOID



Demonstrate confidence.

### SQUARED



Strike an athletic tone.

### TRIANGLE-TAPERED



Exhibit mature experience.